

ASTM D 1177

SCOPE

Standard tests methods for freezing point of engine antifreeze products

FREEZING POINT OF AQUEOUS ANTIFREEZE SOLUTION - REF 9411520 -

MAIN CHARACTERISTICS

- ✓ Freezing point
- ✓ Apparatus of engine antifreeze, complete with thermometer 75F et 76 F

SCOPE OF DELIVERY

FREEZING POINT delivered ready for use with :

REF 23239 Dewar freezing tube (200 ml)

REF 9411522 Dewar cooling bath

REF 15351 ASTM thermometer (75 F), range from -35 to +35°F

REF 15352 ASTM thermometer (76 F), range from -65 to +5°F

NECESSARY ACCESSORY

REF 941597 Dry ice machine



ORDERING INFORMATION

REF 9411520

FREEZING POINT for use on AC 230 V, 50 Hz - 0.3 A

(W) 200x (D) 300x (H) 600 mm (± 6 kg)

CONTACT : sales@normalab.com

NORMALAB FRANCE SAS
ZA Caux Multipôles 1 - F-76 190 Valliquerville
Tel. : +33 232 700 100
Fax : +33 232 704 732

www.normalab.com

REF 9411520 Leaflet 14rev1

DISTRIBUTED BY

G-Labo Germany
Bgm.-Horneffstr.26
69509 Mörlenbach
Tel.: + 49 6209 797100
Fax: + 49 6209 797101
Mail: info@g-labo.de
Web: www.g-labo.de